ENGLISH

Connect this device to the radio transceiver to communicate wirelessly between device and carephone. In order for the device to send alarms through the carephone, you need to pair them.

How to pair the radio transceiver with the carephone

1. Connect the **9350** to the **voice alarm**.
2. Activate the radio learning mode in the Carephone / System
3. **On the voice alarm**; press the red button (Knop logo) to initiate pairing. A confirmation beep is heard on the carephone.
4. **On the carephone**; press the **alarm button** to assign the new device. A new confirmation beep is heard.

For more information on pairing the carephone with the radio transceiver, see the manual of the carephone.

How to test the alarm after pairing

On the voice alarm; press the red button (Knop logo) to generate an alarm. Check that the alarm is received at the Alarm Receiving Center.

This is how it works

1. When an alarm is triggered on the device, the radio transceiver sends the alarm to the carephone.
2. The carephone passes the alarm to the alarm receiving centre or another receiver.

Battery warning

The device automatically alarms at low battery. Replace the battery as soon as possible after low battery level has been indicated.

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