


# User manual

# PU101





## ***Table of contents***

Warnings and safety instructions .....	3
Product description .....	5
Testing and commissioning .....	5
Connection .....	5
Adjustment .....	5
Location .....	6
Service .....	6
Battery alarm .....	6
Changing the battery .....	6
Cleaning .....	7
Spare parts and accessories .....	7
Technical data .....	8
Explanation of symbols and approvals .....	9

## **Warnings and safety instructions**



- The battery in this product must compile with the relevant IEC safety standards for batteries.
- Do not consume or eat the battery as there is a risk of chemical burns.
- Keep new and used batteries away from children.
- If the battery compartment cannot be closed properly, discard the product, and keep it away from children.
- Read the intended use before using.
- The product is not water-resistant unless noted in product specification.

### **Contraindications**

In general, the product cannot be used if:

- The disabled client/patient is mentally ill.
- The disabled client/patient is critically ill.
- The client/patient is unable to use the transmitters.

### **Lifetime after mounting of the device**

The battery, if applies, must be replaced when the battery replacement information appears on the display and according to the user manual.

Lifetime is evaluated in relation to the pressure button. This is made according to the supplier's specification.

### **Lifetime battery (products using battery, only)**

KNOP calculates the life of the batteries. See specification in the user manual.

### **Intended use**

In general, the different variants of the medical device (transmitters/receiver system from KNOP Elektronik) are made as assistive aid for impaired/disabled patients to be able to call for assistance.

The various variants of the medical equipment consist of a combination of products (transmitter and receiver) designed to call for help to impaired/disabled patients; The transmitter-receiver system can be activated in different ways:

- For example, patients who actively do this and are aware that they are calling for help, such as e.g., people with walking difficulties who need help going to the toilet.
- Or disabled patients who are not aware that they are inadvertently putting pressure on the sendere.g. during a seizure.



- Or patients with intellectual and cognitive deficits who are not aware that a receiver is receiving a signal from their transmitter when they leave a house or room.

The system is not designed for critically or mentally ill persons.

### **General product description**

The products manufactured by KNOP elektronik consist of several variants of transmitters and receivers that can be combined with each other. In addition, these products are used in combination with positioning and repeater systems.

These products are medical devices intended to call for assistance and are used for disabled/mobile people, such as those with walking difficulties who need help going to the toilet. The transmitters are activated, for example, by a sound or by pressure (e.g. by pressure with a finger or a breath through the mouth). The recipient is supervised by health care personnel or lay people in private homes.

The system is not designed for critically or mentally ill people.

Part of product	Function in the product system
Transmitter	The transmitter can send the signal obtained from the patient to the receiver monitored by the healthcare personnel or lay person. The transmitter products can be activated by button, sound, blow or movements.
Receiver	The transmitters can be coded into all the receivers and to several receivers at the same time. Some receivers also have a summon button for calling assistance from their coworkers.
Repeatersystem	If it is needed to cover a more comprehensive and larger area a repeater system is used. The repeater system also gives an increased functionality as e.g. that an alarm automatically is received at first at the healthcare person closest to the client.
Positionsystem	If a sender is equipped with a position receiver it can be used in connection with a Position system. Not all product variants include position receivers. Typically, it is seen in connection with nursing homes and security for patients with dementia. The receivers can be portable or stationary.

## **Product description**

The product is designed to send codes to KNOP wireless receivers when the user blows on it.

PU101 can be used where the user is unable to activate an alarm with a pull cord or push button.

PU101 must be used in conjunction with a KNOP Elektronik transmitter or existing call systems.


## **Testing and commissioning**

After installation, it is important that the user tests the product by blowing against the propeller.

Staff shall receive an alarm on their alarm system/receiver. Please refer to the user manual of the transmitter or paging system used.

## **Connection**

PU101 has a 3m cord with 3.5mm jack plug.  
Potential-free contact between frame and tip.  
Max. 36V/100mA.



## **Adjustment**

PU101 has the possibility to set multiple sensitivities on the switch inside the product - see photo below.

PU101 can be disassembled without the use of tools. Use a small screwdriver to adjust the switch.



- 1 = High sensitivity:  
Short weak puff.
- 2 = Medium sensitivity  
Slightly longer and stronger puff.
- 3 = Low sensitivity.  
Long strong puff.

0, 4-9: Reserved for other functions.

## Location

PU101 should be positioned so that the user can easily blow straight into the propeller.

The push contact can be mounted on e.g. a gooseneck, contact your retailer for more information.

## Service


When the user blows against the propeller, strongly enough (depending on the sensitivity described earlier) to activate PU101, it will flash green.

## Battery alarm

In case of low battery, a battery alarm with a red light appears through the propeller and on the sides when PU101 is activated with a puff.

## Changing the battery

The top and bottom of PU101 can be unscrewed without tools.



The battery is gently pressed out of the holder, for example with a cotton swab or a match.

**Never use a metal object.**

Carefully insert the new battery into the holder **WITHOUT** using tools, as this may short-circuit the battery.

**Use a CR1632 battery.**

**IMPORTANT + must be facing up.**

After changing the battery, check PU901 as described under "Testing and commissioning". It takes about 1 minute after battery change before it can be activated.


## ***Cleaning***

The product can be cleaned with a moist cloth or disinfected wet wipe.

The propeller can be gently pulled out with pliers or other tools and cleaned.

The product is then thoroughly tested.

If necessary, the propeller can be replaced with a new one. See section "Spare parts and accessories".



## ***Spare parts and accessories***

Spare parts and accessories can be ordered on our website [www.knop.dk](http://www.knop.dk)



## Technical data

<b>Battery type:</b>	1 x CR1632, Lithium 3V 140mAh.
<b>Battery life (expected):</b>	Up to 1 year at 10 activations per day <sup>(1)</sup>
<b>Low battery alarm:</b>	At approx. 1/3 remaining capacity
<b>Power consumption (inactive):</b>	<1µA.
<b>Power consumption (active):</b>	<7mA for 3 seconds per activation.
<b>Ambient environment:</b>	Indoor use ≤ 90% <b>non-condensing</b>
<b>Ambient temperature:</b>	0°C to +40°C
<b>Enclosure type:</b>	White polyoxymethylene
<b>Enclosure dimensions:</b>	Ø: 43mm L: 106mm
<b>Density:</b>	IP40
<b>Weight:</b>	115g

The right to make changes is reserved.















All rights reserved.

**KNOP** Elektronik A/S

<sup>(1)</sup> Expected battery life depends on usage and battery quality.



## Explanation of symbols and approvals

	<p><b>This product complies with:</b></p> <table border="0"> <tr> <td>Directive 2017/745/EU</td><td>MDR</td></tr> <tr> <td>Directive 1907/2006/EU</td><td>REACH</td></tr> <tr> <td>Directive 2011/65/EU</td><td>RoHS</td></tr> <tr> <td>Directive 2012/19/EU</td><td>WEEE</td></tr> <tr> <td>ISO 14971:2019</td><td>Risk Management for Medical devices</td></tr> <tr> <td>EN 301 498-1 V2.2.3</td><td>Electro Magnetic Compatibility</td></tr> <tr> <td>EN 301 489-3 V3.1.1</td><td>Electro Magnetic Compatibility</td></tr> <tr> <td>EN 50130-4:2011 + 2014</td><td>Immunity alarm systems</td></tr> <tr> <td>EN 300 220-1 V3.1.1</td><td>Short Range Devices</td></tr> <tr> <td>EN 300 220-2 V3.1.1</td><td>Short Range Devices</td></tr> <tr> <td>EN 300 220-3 V2.1.1</td><td>Short Range Devices</td></tr> <tr> <td>EN 60601-1-2:2014 + 2015</td><td>Electromagnetic compatibility (Medical)</td></tr> <tr> <td>EN 62368-1:2020</td><td>Electrical safety</td></tr> <tr> <td>EN 50581:2012</td><td>Hazardous substances</td></tr> </table>	Directive 2017/745/EU	MDR	Directive 1907/2006/EU	REACH	Directive 2011/65/EU	RoHS	Directive 2012/19/EU	WEEE	ISO 14971:2019	Risk Management for Medical devices	EN 301 498-1 V2.2.3	Electro Magnetic Compatibility	EN 301 489-3 V3.1.1	Electro Magnetic Compatibility	EN 50130-4:2011 + 2014	Immunity alarm systems	EN 300 220-1 V3.1.1	Short Range Devices	EN 300 220-2 V3.1.1	Short Range Devices	EN 300 220-3 V2.1.1	Short Range Devices	EN 60601-1-2:2014 + 2015	Electromagnetic compatibility (Medical)	EN 62368-1:2020	Electrical safety	EN 50581:2012	Hazardous substances
Directive 2017/745/EU	MDR																												
Directive 1907/2006/EU	REACH																												
Directive 2011/65/EU	RoHS																												
Directive 2012/19/EU	WEEE																												
ISO 14971:2019	Risk Management for Medical devices																												
EN 301 498-1 V2.2.3	Electro Magnetic Compatibility																												
EN 301 489-3 V3.1.1	Electro Magnetic Compatibility																												
EN 50130-4:2011 + 2014	Immunity alarm systems																												
EN 300 220-1 V3.1.1	Short Range Devices																												
EN 300 220-2 V3.1.1	Short Range Devices																												
EN 300 220-3 V2.1.1	Short Range Devices																												
EN 60601-1-2:2014 + 2015	Electromagnetic compatibility (Medical)																												
EN 62368-1:2020	Electrical safety																												
EN 50581:2012	Hazardous substances																												
	<p>The product must not be used if the packaging is damaged.</p>																												
	<p><b>Medical Device</b> Class 1, rule 1</p>																												
	<p><b>Manufacturer</b> KNOP Elektronik A/S, Fabriksvej 20, DK-7600 Struer, Denmark</p>																												
	<p>Read the manual(s) before installation and commissioning at <a href="http://www.knop.dk">www.knop.dk</a>.</p>																												
	<p>Interference may occur in an environment with equipment marked with this symbol.</p>																												
	<p>Must be protected against liquids.</p>																												
	<p>0 °C to +40 °C, temperature limit for transport/storage and use.</p>																												
	<p>The product must not be disposed of with normal household waste.</p>																												
	<p><b>Single Registration Number</b> DK-MF-000025631</p>																												
	<p><b>Unique Device Identifier</b> PU101: 05744002852032</p>																												
	<p><b>Product reference</b> PU101: Blow switch</p>																												
	<p><b>Serial number</b> Placed on the product</p>																												
	<p><b>Ingress Protection code</b> IP40</p>																												