# Extreme temp 9350

## ENGLISH

In this device we have integrated a radio transceiver to communicate wirelessly between device and carephone. In order for the device to send alarms through the carephone, you need to pair them.

### Preparations

1. Disassemble the device.
2. Insert the included CR2032 coin cell battery with "+" pointing up.
3. Reassemble the device.

### How to pair the radio transceiver with the carephone

1. **On the carephone**; press and hold the **green button** for four seconds until a beep sounds to start pairing mode. Release the button.
2. **On the Extreme temp**, press the test button to generate an alarm. A confirmation beep is heard on the carephone.
3. **On the carephone**; press the **alarm button** to assign the new device. A new confirmation beep is heard.

For more information on pairing the carephone with the radio transceiver, refer to the manual of the carephone.

### How to test the alarm after pairing

**On the Extreme temp**; press the test button to generate an alarm. Check that the alarm is received at the Alarm Receiving Center.

*Note! The radio transceiver is blocked for further alarms for 10 seconds after an alarm is generated.*

### This is how it works

1. When an alarm is triggered on the device, the radio transceiver sends the alarm to the carephone.
2. The carephone passes the alarm to the alarm receiving centre or another receiver.

### Battery warning

The device automatically alarms at low battery. Replace the battery as soon as possible after low battery level has been indicated.

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©2025 Careium Sweden AB. All rights reserved.